



Impression Transfer

Available for both open-tray and closed-tray technique.




Ball Attachments

Conical and Hex Connection manufactured in Titanium Alloy.




Healing Abutment

Available in various heights and diameters.



Implant Prosthetic Indicator (IPI)

Conical and Hex Connection manufactured in Titanium Alloy.




Cement-Retained Abutments

Anatomic and Angulated Zirconia Abutments with Conical or Hex Connections.

Anatomic Zirconia with Titanium Base
Base Ø4.5mm with Length 1mm, 2mm, 3mm

15° Angulated Zirconia with Titanium Base
Base Ø5.35mm with Length 1mm, 2mm, 3mm



www.SpiralTech.com

OUR COMPANY



Founded in 2012, SpiralTech was created to fill a void in the dental implant world, simplified treatment and paired with the best technologies and service.

We wanted to provide innovative designs, along with the necessary education and service, to doctors of all experience levels.

Ultimately, we desired to enhance the doctor and patient experience. The ESi implant product line, is the first step in what we see as a continuous journey towards clinical excellence.



Prosthetics & Abutments

HEADQUARTERS



875 N. Michigan Avenue
Chicago, IL 60611
1 (312) 440-7777
Toll free: 1 (888) 991-7997

ISO 13485:2003
EC Certificate US15/842415
F/7.5.1.1_49 Rev.03 3.3.17

© 2016 SpiralTech Superior Dental Implants, Inc. All Rights Reserved.



Cement-Retained Abutments

Prefabricated Flat & Anatomic Abutments with Conical or Hex Connections.


Hex (for Ø3.5, Ø4.3, Ø5.0, Ø6.0 implants)

NP (for Ø3.5 implants)


RP (for Ø4.3, Ø5.0, Ø6.0 implants) with shoulder/cuff collar design

Design	Diameter (mm)	Facial Height (mm)	Hex Catalog No.	NP Conical Catalog No.	RP Conical Catalog No.	
Flat, titanium, Narrow	3.8	5	FTAH3805	FTACNP3805	FTACRP3805	
		7	FTAH3807	FTACNP3807	FTACRP3807	
		9	FTAH3809	FTACNP3809	FTACRP3809	
		11	FTAH3811	FTACNP3811	FTACRP3811	
Flat, titanium, Standard	4.5	5	FTAH4505	FTACNP4505	FTACRP4505	
		7	FTAH4507	FTACNP4507	FTACRP4507	
		9	FTAH4509	FTACNP4509	FTACRP4509	
		11	FTAH4511	FTACNP4511	FTACRP4511	
Flat, titanium, Wide	5.5	5	FTSAH5505		FTSARP5505	
		7	FTSAH5507		FTSARP5507	
		9	FTSAH5509		FTSARP5509	
Flat, titanium, with Straight Shoulder	3.8	1	FTASH3810	FTASCNP3810	FTSACRP3810	
		2	FTASH3820	FTASCNP3820	FTSACRP3820	
		3	FTASH3830	FTASCNP3830	FTSACRP3830	
	4.5	1	FTASH4510	FTASCNP4510	FTSACRP4510	
		2	FTASH4520	FTASCNP4520	FTSACRP4520	
		3	FTASH4530	FTASCNP4530	FTSACRP4530	
		4	FTASH4540	FTASCNP4540	FTSACRP4540	
	5.5	1	FTASH5510		FTSACRP5510	
		2	FTASH5520		FTSACRP5520	
		3	FTASH5530		FTSACRP5530	
	Anatomic, titanium	4.5	1	ATAH4510	ATACNP4510	ATACRP4510
			2	ATAH4520	ATACNP4520	ATACRP4520
3			ATAH4530	ATACNP4530	ATACRP4530	
4			ATAH4540	ATACNP4540	ATACRP4540	
5.5		1	ATAH5510		ATACRP5510	
		2	ATAH5520		ATACRP5520	
		3	ATAH5530		ATACRP5530	
		Angulated, titanium 15°	4.7	1	15ATAH4710	15ATACNP4710
2	15ATAH4720			15ATACNP4720	15ATACRP4720	
3	15ATAH4730			15ATACNP4730	15ATACRP4730	
Angulated, titanium 25°	5.2	1	25ATAH5210	25ATACNP5210	25ATACRP5210	
		2	25ATAH5220	25ATACNP5220	25ATACRP5220	
		3	25ATAH5230	25ATACNP5230	25ATACRP5230	


Design	Diameter (mm)	Facial Height (mm)	Slim Catalog No.
Flat titanium, slim	3.7	1	15ASACNP3710
	3.7	2	15ASACNP3720
	3.7	3	15ASACNP3730



Flat, Narrow, Titanium



Flat, Standard, Wide, Titanium



Flat, Titanium with Shoulder



Anatomic, Titanium



15° Angulated Titanium



25° Angulated Titanium